

TOMATO GLASSES



EYEWEAR OF YEAR 2016 AWARD WINNER (JAPAN)
 OPTRA FAIR AWARD WINNER 2013 (UK)
 KIDS DESIGN AWARD 2014 (JAPAN)
 GOOD DESIGN (JAPAN)
 GOOD DESIGN (KOREA)
 ODMA AWARDS 2017 (AUSTRALIA)
 O=MEGA AWARDS 2019 WINNER (AUSTRALIA)
 ONE MILLION DOLLAR, PRIZE FOR TOP OF EXPORT
 THREE MILLION DOLLAR, PRIZE FOR TOP OF EXPORT

I wear you ♥ You frame me ♥ I love you

OUR STORY

Our journey started when I noticed that my 24 month old son's eyes were not aligning correctly. Our trip to the ophthalmologist confirmed Accommodative Esotropia.

His condition Hyperopia is a result of a small eye, which caused his eyes to strain so much that they crossed and eventually one eye became lazy.

The only correction was to wear thick convex lenses and to wear an eye patch to prevent a lazy eye known as Amblyopia.

The limited choice of frames on the market at the time with nose pads were only in metal.

They were heavy and bulky and constantly kept falling off his face. Compliance was poor as my son refused to wear them due to the discomfort and weight. His condition worsened.

There was nothing available that could cater for a small child's head with a flat bridge. My disappointment was such that I was compelled to make his frames on my own.

This sparked the beginning of Tomato Glasses : a company and product born out of the affection and love I have for my child. His vision improved and with the help of surgery his vision is now normal.

I am glad I did not give up on him and ever so grateful that the Tomato Glasses helped him. I am also humbled to know that my passion and love for my son is now helping and supporting other children and parents across the globe that require custom fitted, lightweight and comfortable glasses. Children can now leave these spectacles on, with amazing success and great comfort. Tomato Glasses are frames that fit to a child's face rather than the child fitting the frame. We, Tomato Glasses, will continue to make children's frames so that we can make children and their parents happy.



Scan For Video



NEW MASCOT

Non-Slip! Safe! Adjustable!

TOMATO GLASSES FEATURES

01 DURABLE GLASSES

1. The frames are made using non-toxic materials that are approved by the FDA.
2. The frames are very flexible, display excellent elasticity and do not break easily.



02 COMFORTABLY FITTING GLASSES

1. The frames are incredibly light: only 6 to 11g in weight.
2. There is a cushion on the tip of the frames that matches the complex curves of the ears and will not create pressure points.



Weight by type

Baby 6.3g | Kids 7.0-9.5g | Junior 8.5-10g | Young 9.4-11.5g

03 NON SLIP GLASSES

The bridge of the glasses matches the curve of the nose and therefore holds the glasses comfortably in place and prevents them from slipping down.




04 ADJUSTABLE TEMPLE LENGTHS

It is easy and simple to adjust the length of the temples. To shorten them, remove the ear tip and cut off the excess length. Replace ear tips and secure these in position using the screws. To lengthen the temples, pull the ear tips back and secure in place with the screws.



05 ADJUSTABLE NOSE PAD HEIGHTS

The nose pads can be attached to three points of the frame (five points on the baby range frames). As every child's face is unique, this feature allows the correct placement and height to be achieved for each child.



06 A STANDARD HEAD STRAP

All Tomato Glasses frames come with a head strap which attaches to the ear tips. They have been designed to prevent the glasses from slipping even with high powered lenses. The head strap provides security when running and jumping, making Tomato Glasses great for active children. They are also perfect for children with sensitive ears.



Baby A



A very flexible frame for the utmost comfort which is highly adjustable to accommodate a baby's rapid growth.



SIZE 35□13 37□13 39□14

FRAME SIZE	35	37	39
NOSE PAD TYPE	B1	B1	B1

TBAC1



TBAC2



TBAC3




SPECIAL FEATURES

POINT 1



There are 5 holes on the front frame to attach the nose pad.

POINT 2



There are 3 different sizes of nose pads that can decide the distance between the eyes and lenses.

POINT 3



Safe material, flexible and adjustable temples.

Kids A

OVAL FRAMES



SIZE 37□14 40□14 43□15 45□16

FRAME SIZE	37	40	43	45
NOSE PAD TYPE	A1	A1	W1	W2

TKAC1



TKAC13



TKAC3



TKAC14



TKAC4



TKAC15



TKAC7



TKAC25



TKAC11



TKAC26



KIDS A

SIZE 37□14 40□14 43□15 45□16

TKAC27



TKAC32



TKAC28



TKAC33



TKAC29




TKAC34(Matte)



SIZE 37□14 40□14 43□15 45□16

FRAME SIZE	37	40	43	45
NOSE PAD TYPE	A1	A1	W1	W2

TKAC36



blue glitter




NEW

TKAC39



TKAC37




Same flexible material as Baby line



TKAC509



TKAC38



Same flexible material as Baby line



Kids B

CIRCULAR FRAMES



SIZE 40□14 42□15 44□16 46□17

FRAME SIZE	40	42	44	46
NOSE PAD TYPE	A1	W1	W1	W2

TKBC1



TKBC7



TKBC2



TKBC8



TKBC3



TKBC9



TKBC4



TKBC11



TKBC5



TKBC14




KIDS B

SIZE 40□14 42□15 44□16 46□17

TKBC18



TKBC26



blue glitter

TKBC21

(Dark navy)



NEW

TKBC27(Mat)



TKBC22

(navy)



NEW

TKBC28



TKBC23



NEW

TKBC501(Glossy)



TKBC24(Mat)



Kids C

RECTANGULAR FRAMES



KIDS C

SIZE 40□14 42□15 44□16 46□17

TKCC1



TKCC6

(Dark navy)



TKCC2



TKCC7



TKCC3



TKCC8



TKCC4

(Crystal grey)



TKCC9



TKCC5



TKCC14

(Black)



SIZE 40□14 42□15 44□16 46□17

FRAME SIZE	40	42	44	46
NOSE PAD TYPE	A1	W1	W1	W2

TKCC16
(Navy)



TKCC18




TKCC21



TKCC22(Matte)



TKCC24



blue glitter

NEW

TKCC25(Matte)



NEW

TKCC26



NEW

TKCC501(Glossy)



Kids D

RECTANGULAR FRAMES

 **NEW SIZE 46 has been added.**



NEW

SIZE 38□17 40□17 42□17 44□17 46□17

FRAME SIZE	38	40	42	44	46
NOSE PAD TYPE	A1	A1	W1	W1	W2

TKDC1



TKDC8

(Dark navy)



TKDC2



TKDC11



TKDC3



TKDC13



TKDC5



TKDC14



TKDC7

(Black)



TKDC17

(Crystal brown)



KIDS D

SIZE 38□17 40□17 42□17 44□17 **46□17**

NEW

TKDC18



TKDC24



TKDC19



TKDC25



TKDC21 (Matte)



NEW

TKDC501 (Glossy)



TKDC22



Junior A

RECTANGULAR FRAMES



JUNIOR A

SIZE 47□17 49□17 51□17

FRAME SIZE	47	49	51
NOSE PAD TYPE	W2	W3	W3

TJAC1



TJAC11



TJAC2



TJAC12



TJAC3



TJAC13



TJAC4



TJAC16(Matte)




TJAC5

(Crystal coffee)



TJAC17



blue glitter

TJAC6



Junior B

Discontinued

As of November 2022, TJBC(Junior B) line has been discontinued.

You are able to order only parts of TJBC from our stock when needed.

Our new style Junior D has been released as a replacement of Junior B.



TJBC4



TJBC12



TJBC13



TJBC5



TJBC14



TJBC6



TJBC15



TJBC7



TJBC16



TJBC9



TJBC17(Matte)



Junior C

RECTANGULAR FRAMES



7yrs to
adults



SIZE 46□17 48□17 50□17 52□17 54□17.5

FRAME SIZE	46	48	50	52	54
NOSE PAD TYPE	W2	W2	W3	W3	W3

TJCC1



TJCC8

(Black)



TJCC2



TJCC11



TJCC3




TJCC13(Matte)



TJCC4



TJCC14



blue glitter

TJCC6



NEW

TJCC15



TJCC7



Junior D



7yrs to
adults

A round shape frame which looks great on anyone
Suitable frame for Myopia Progression Control Lenses



SIZE 46□17 48□17 50□18 52□18

FRAME SIZE	46	48	50	52
NOSE PAD TYPE	W2	W2	W3	W3

TJDC2



TJDC9(Glossy)



TJDC3



TJDC11(Matte)



TJDC4



NEW

TJDC12



TJDC5



NEW

TJDC13



TJDC7



NEW

TJDC14



TJDC8



NEW

TJDC15



Young A


For late teens to adults

Frame for High Prescription Lenses

Resolve the appearance problems caused by the high prescription lenses!

Resolve eyeglasses' thickness and weight problems caused by the high prescription lenses!


Keep all the features of Tomato Glasses, especially about the non-slipping characteristic!



Young A frame for High Prescription Lenses



Resolve the appearance problems caused by the high prescription lenses!



Resolve eyeglasses' thickness and weight problems caused by the high prescription lenses!



Keep all the features of Tomato Glasses, especially about the non-slipping characteristic!



SIZE S(37□27) M(39□29) L(41□31)

FRAME SIZE	S(37□27)	M(39□29)	L(41□31)
NOSE PAD TYPE	W2	W2	W3

TYAC1



TYAC4



TYAC2



TYAC5



TYAC3



TYAC6(Matte)



Weight felt on top of nose




-5.00D Lenses



-10.00D Lenses

Frame Sizes



Recommended Frame Sizes

PD(mm)	Frame Size(FPD)
59, 60, 61, 62	S (64)
63, 64, 65, 66	M (68)
67, 68, 69, 70	L (72)



Young B



For late teens
to adults

Chic and sophisticated design
A double frame with high durability and refined texture
Suitable for wide face (size 56mm frames)



SIZE 50□19 53□19 56□19

FRAME SIZE	50	53	56
NOSE PAD TYPE	W3	W3	W3

TYBC1



TYBC5



TYBC2



TYBC6(Matte)



TYBC3



TYBC7



TYBC4



Young B special case includes both folding and pouch cases. One of each A4 and W3 nose pads are included for extra.



NOSE PAD FITTING KIT


Further research has indicated the need to cater for a variety of nose shapes. Therefore, we have developed a range of different sized nose pads, Baby fit (B1, B2, B3), Asian fit (A1, A2, A3, A4, A5) and Wide fit (W1, W2, W3).


All nose pads are interchangeable on frames and the fitting kit stores all nine sizes of nose pad.

This kit will allow optometrists and opticians to custom fit the frames to each child, therefore, offering the best customer service.




Nose Pad Fitting Kit

B SERIES (BABY)	 <p>B1 TYPE</p> <ul style="list-style-type: none"> ① external width 22mm ② internal width 8mm ③ height 8mm  <p>5mm</p>	A SERIES	 <p>A1 TYPE</p> <ul style="list-style-type: none"> ① external width 25mm ② internal width 10mm ③ height 7mm 	W SERIES	 <p>W1 TYPE</p> <ul style="list-style-type: none"> ① external width 23mm ② internal width 11mm ③ height 8mm
	 <p>B2 TYPE</p> <ul style="list-style-type: none"> ① external width 22mm ② internal width 8mm ③ height 9mm  <p>6.5mm</p>		 <p>A2 TYPE</p> <ul style="list-style-type: none"> ① external width 25mm ② internal width 10mm ③ height 7mm 		 <p>W2 TYPE</p> <ul style="list-style-type: none"> ① external width 25mm ② internal width 13mm ③ height 8mm
	 <p>B3 TYPE</p> <ul style="list-style-type: none"> ① external width 22mm ② internal width 8mm ③ height 10mm  <p>8mm</p>		 <p>A3 TYPE</p> <ul style="list-style-type: none"> ① external width 27mm ② internal width 12mm ③ height 7mm 		 <p>W3 TYPE</p> <ul style="list-style-type: none"> ① external width 26mm ② internal width 15mm ③ height 8mm
		A SERIES	 <p>A4 TYPE</p> <ul style="list-style-type: none"> ① external width 25mm ② internal width 10mm ③ height 7mm 		
			 <p>A5 TYPE</p> <ul style="list-style-type: none"> ① external width 25mm ② internal width 10mm ③ height 7mm 		




THE RELATIONSHIP BETWEEN THE TYPES OF NOSE PADS AND FITTING RESULTS



If the frame sits too far away from a child's eyes, a wider nose pad should be used (e.g. replace a W1 for a W3 nose pad).
 If the frame sits too close to eyelashes then the above spacer kit can be used to create a greater space.


FOR INFANTS (B1,B2,B3), PLEASE REPLACE IT WITH THE LOWER NUMBER.

FUNCTION OF A4, A5 NOSE PADS



The A4 and A5 nose pads sets raise the eyeglass frame further up.

SPACER KIT



1 2mm type	2 3mm type	3 4mm type

Our Research and Development team have created a Spacer Kit.

It is designed to create space between the lenses and a child's eyelashes or face when using high powered lenses.

CONTENTS

- 1 2mm / 2 set
- 2 3mm / 2 set
- 3 4mm / 2 set



TOMATO GLASSES CASSETTE AND ACCESSORIES



- Accessories are included in the Tomato Glasses cassette.
- The Baby frames have more accessories in their packs.
- Very convenient for storing parts.

Contents

- case · frame · strap · extra nose pads
- extra tips · extra tip screws
- driver · fabric cleaner · user guide

EYE CATCHING PROMOTIONAL ITEMS

※ Contents and design are subject to change without any notice.



Vinyl banner size 60x160cm



Poster size 42.6x60.6cm



Acrylic tabletop sign size 10x14cm



Soft case



Foam board size 36x25cm



Acrylic tabletop sign size 35.5x26.5cm

Wood Display Tray



1set : 4ea



Plastic Display Tray size 49.5x41cm

OUR AWARDS

TOMATO GLASSES are an award winning specialized brand of children's glasses. We supply frames worldwide and we have a good reputation for design and quality from ophthalmologists and opticians around the world.

WINNER



AWARDS OF EXCELLENCE

O-MEGA AWARDS 2019 WINNER (AUSTRALIA)



GOOD DESIGN AWARD 2017

GOOD DESIGN (JAPAN)

FINALIST



AWARDS OF EXCELLENCE

ODMA AWARDS 2017 (AUSTRALIA)



EYEWEAR OF THE YEAR 2016 AWARD WINNER (JAPAN)



KIDS DESIGN AWARD 2014 (JAPAN)



GOOD DESIGN
GOOD DESIGN (KOREA)



OPTRAFAIR
Birmingham 2013
WINNER
OPTRA FAIR AWARD WINNER 2013(UK)



One million dollar, Prize for Top of export



Three million dollar, Prize for Top of export



EYEWEAR OF TEAR 2016 AWARD WINNER (JAPAN)



KIDS DESIGN AWARD 2014



Optra Awards 2013 1st



Certificate of Good Design Selection



ODMA AWARDS 2017 (AUSTRALIA)



GOOD DESIGN (JAPAN)



AUSTRALIA O-MEGA AWARD



Toxin-Not Detected TEST REPORT(ear tips)



Toxin-Not Detected TEST REPORT(nose pads)



TEST REPORT(frame) Toxin-Not Detected



Certificate Europe



No.10-1179157 Certificate of patent



No.10-1186259 Certificate of patent



Certificate of Trade Service Mark Registration



Certificate of Trademark Registration (JAPAN)



Certificate of Trademark Registration (USA)



No.30-0652187 Certificate of Design



No.30-0652189 Certificate of Design



No.30-0652192 Certificate of Design



Certificate of Origin




No.30-1027921 Certificate of Design



Certificate of Trademark Registration (INDIA)



Certificate of Trademark Registration (Europe)



HEAD OFFICE

HEAD OFFICE (Korea) : TOMATO GLASSES (+82) 51-244-8760

DISTRIBUTORS

- Australia : TOMATO GLASSES AUSTRALIA (+61) 3-9448-8932
- Brazil : Vis-à-vis Importação (+55) 11-26797902
- Bulgaria : Milde Eood (+359) 888-811-508
- Canada : TOMATO GLASSES CANADA (+1) 778-954-7880
- China : Shanghai tomato glassesCo.,Ltd (+86) 189-1281-2032
JiangSu Tomato Glasses Co.,Ltd (+86) 139-1395-7255
GuangZhou Xin die cui glasses com.,Ltd (+86) 139-2217-6767
- Colombia : TIENDA OPTICA (+57) 909-485-5592
- Czech Republic : TOMATO GLASSES Czech Republic (+420) 605747864
- Greece : F-ORASIS (+30) 2310-445-500
- HongKong : WING KI ENTERPRISE LIMITED (+852) 6014-6863
- Indonesia : PT.Optik Klinik Mata Nusantara (+62) 21-7581-6601
- Japan : TOMATO GLASSES 092-263-8166
- Malaysia : TOMATO GLASSES (SINGAPORE) PTE LTD (+65) 9022-8960
- Middle East : TOMATO GLASSES ME (+971) 4-252-4596
- New Zealand : TOMATO GLASSES NZ (+64) 21-869-300
- Philippines : TOMATO GLASSES PH (+63) 920-932-6578
- Singapore : TOMATO GLASSES (SINGAPORE) PTE LTD (+65) 9022-8960
- South Africa : Danker Lenses (+27) 11-262-3980
- Taiwan : ABBY Optical Limited (+886) 4 24711616 / GinDin international Co.,Ltd (+886) 2-2273-7785
- Türkiye : MK OPTİK İÇ VE DIŞ TİCARET LİMİTED ŞİRKETİ (+90) 312-467-27-12
- United Kingdom : TOMATO GLASSES UK / TOMATO GLASSES EU (+44) 1395-488188
- United States : TOMATO GLASSES USA (+1) 201-655-5758

All employees at Tomato Glasses are committed to better vision for infants and children.


All of our employees are diligently working to deliver new ideas to improve vision for infants and children.


“Let every child see the world clearly.”

“Every child should be able to study, to exercise, and to play freely without any restraints.”

“Every child is precious, and we hope every child to grow up healthy.”

With these kinds of thoughts in mind, our employees are relentlessly striving together to make better products.





TOMATO GLASSES

HEAD OFFICE

29 Gudeok-ro 120 beon-gil, Seo-gu, Busan 49246, Korea

Tel : +82-51-244-8760, Fax : +82-51-244-8761

www.tomatoglasses.com

FACTORY

21, Haeyang-ro 57beon-gil, Youngdo-gu, Busan 49007, Republic of Korea

Tel : +82-70-7776-8760

JAPAN OFFICE

2-20 Hakatachuo Bld.2F, Nakagofukumachi, Hakata-ku,

Fukuoka-shi, Fukuoka, 812-0035, Japan

Tel : +81-92-263-8166, Fax : +81-92-263-8167

www.tomatoglasses.co.jp

