







I wear you You Frame me I love you









rmingham 20° **WINNER** OPTRA FAIR AWARD WINNER 2013 (UK)



KIDS DESIGN AWARD 2014



















OUR STORY

Our journey started when I noticed that my 24 month old son's eyes were not aligning correctly. Our trip to the ophthalmologist confirmed Accommodative Esotropia.

His condition Hyperopia is a result of a small eye, which caused his eyes to strain so much that they crossed and eventually one eye became lazy.

The only correction was to wear thick convex lenses and to wear an eye patch to prevent a lazy eye known as Amblyopia.

The limited choice of frames on the market at the time with nose pads were only in metal.

They were heavy and bulky and constantly kept falling off his face. Compliance was poor as my son refused to wear them due to the discomfor \Box t and weight. His condition worsened.

There was nothing available that could cater for a small child's head with a flat bridge.

My disappointment was such that I was compelled to make his frames on my own.

This sparked the beginning of Tomato Glasses: a company and product born out of the affection and love I have for my child. His vision improved and with the help of surgery his vision is now normal.

I am glad I did not give up on him and ever so grateful that the Tomato Glasses helped him.

I am also humbled to know that my passion and love for my son is now helping and supporting

other children and parents across the globe that require custom fitted, lightweight and comfortable glasses. Children can now leave these spectacles on, with amazing success and great comfort. Tomato Glasses are frames that fit to a child's face rather than

the child fitting the frame.

We, Tomato Glasses, will continue to make children's frames so that we can make children and their parents happy.



Scan For Video









Durability Test

Bridge deformation Test

High temperature stability Test







ISO 9001 Certification

ISO 13485 Certification




Non-Slip! Safe! Adjustable!



TOMATO GLASSES FEATURES

11 DURABLE GLASSES

- 1. The frames are made using non-toxic materials that are approved by the FDA.
- 2. The frames are very flexible, display excellent elasticity and do not break easily.







(1) COMFORTABLY FITTING GLASSES

- 1. The frames are incredibly light: only 6 to 11g in weight.
- 2. There is a cushion on the tip of the frames that matches the complex curves of the ears and will not create pressure points.





Weight by type

Baby 6.3-6.8g | Kids 7.0-9.5g | Junior 8.5-10g | Young 9.4-11.5g | Active 8.5-10g

03 NON SLIP GLASSES

The bridge of the glasses matches the curve of the nose and therefore holds the glasses comfortably in place and prevents them from slipping down.



14 ADJUSTABLE TEMPLE LENGTHS

It is easy and simple to adjust the length of the temples.

To shorten them, remove the ear tip and cut off the excess length. Replace ear tips and secure these in position using the screws. To lengthen the temples, pull the ear tips back and secure in place with the screws.



05 ADJUSTABLE NOSE PAD HEIGHTS

The nose pads can be attached to three points of the frame (five points on the baby range frames).

As every child's face is unique, this feature allows the correct placement and height to be achieved for each child.





06 A STANDARD HEAD STRAP

All Tomato Glasses frames come with a head strap which attaches to the ear tips. They have been designed to prevent the glasses from slipping even with high powered lenses.

The head strap provides security when running and jumping, making Tomato Glasses great for active children. They are also perfect for children with sensitive ears.





which is highly adjustable to accommodate a baby's rapid growth.









 FRAME SIZE
 35
 37
 39

 NOSE PAD TYPE
 B1
 B1
 B1




SPECIAL FEATURES

POINT 1



There are 5 holes on the front frame to attach the nose pad.

POINT 2



There are 3 different sizes of nose pads that can decide the distance between the eyes and lenses.

POINT 3



Safe material, flexible and adjustable temples.



The Kids A model has been discontinued as of November 2025

However, replacement parts are still in stock and can be ordered as needed

- * A replacement model, Kids E, is available
- * Kids A (TKAC) and Kids E (TKEC) are not interchangeable
- * The actual product color may differ from the image due to printing































Kids B

CIRCULAR FRAMES

* Use the same temple as Kids C (TKCC)










 FRAME SIZE
 40
 42
 44
 46


 NOSE PAD TYPE
 A1
 W1
 W1
 W2



KIDS B

SIZE 40\(\)14\(\)42\(\)15\(\)44\(\)16\(\)46\(\)17\(\)

FRAME SIZE 40 42 44 46 W1 W1 W2 NOSE PAD TYPE A1









Kids C

RECTANGULAR FRAMES

* Use the same temple as Kids B (TKBC)










KIDS C SIZE 40014 42015 44016 46017











Kids D

RECTANGULAR FRAMES

46 has been added.


* Use the same temple as Kids F(TKFC)












KIDS D SIZE 38-17 40-17 42-17 44-17 46-17











Kids E

A simple oval-shaped model designed with images children love

* Although it looks similar to the Kids A (TKAC),
its frame front and temples are not compatible
with Kids A















| FRAME SIZE | 37 | 40 | 43 | 45 | 47 |
|---------------|----|----|----|----|----|
| NOSE PAD TYPE | A1 | A1 | W1 | W2 | W2 |









Kids F


Featuring a soft, rounded frame with temples decorated in various designs, this model is also ideal for functional eyeglass lenses designed to help prevent myopia

* Use the same temple as Kids D(TKDC)



























Teenagers and older

Junior C

Simple square design with an intellectual feel

















Teenagers and older

Junior D

A round shape frame which looks great on anyone
Suitable frame for Myopia Progression

Control Lenses
















Young A



Resolve the appearance problems caused by the high prescription lenses!

Resolve eyeglasses' thickness and weight problemscaused by

the high prescription lenses!

Keep all the features of Tomato Glasses, especially about the non-slipping characteristic!









Young A frame for High Prescription Lenses










Resolve the appearance problems caused by the high prescription lenses!

Resolve eyeglasses' thickness and weight problems caused by the high prescription lenses! Keep all the features of Tomato Glasses, especially about the non-slipping characteristic!

 FRAME SIZE
 S(37□27)
 M(39□29)
 L(41□31)

 NOSE PAD TYPE
 W2
 W2
 W3




Weight felt on top of nose





-10.00D Lenses

Frame Sizes



Recommended Frame Sizes

| PD(mm) | Frame Size(FPD) | | |
|----------------|-----------------|--|--|
| 59, 60, 61, 62 | S (64) | | |
| 63, 64, 65, 66 | M (68) | | |
| 67, 68, 69, 70 | L (72) | | |









Teenagers and older

Young B

Chic and sophisticated design

A double frame with high durability and refined texture

Suitable for wide face (size 56mm frames)

The Young series includes an extra nose pad for each type, A4 and W3.









| FRAME SIZE | 50 | 53 | 56 |
|---------------|----|----|----|
| NOSE PAD TYPE | W3 | W3 | W3 |









Teenagers and older

Young C

Unique double-frame for a distinctive look

56mm option for larger faces

The Young series includes an extra nose pad for each type, A4 and W3.










| FRAME SIZE | 50 | 53 | 56 |
|---------------|----|----|----|
| NOSE PAD TYPE | W3 | W3 | W3 |









Teenagers and older

Active A

For those who prefer a simple square design

Recommended for comfort seekers

Ideal for functional lenses









SIZE 50□17 52□17 54□17.5

 FRAME SIZE
 50
 52
 54

 NOSE PAD TYPE
 W3
 W3
 W3









Teenagers and older

Active B

Oval shape with a soft feel Lightweight and comfortable fit Ideal for functional lenses










| FRAME SIZE | 50 | 52 | 54 |
|---------------|----|----|----|
| NOSE PAD TYPE | W3 | W3 | W3 |




NOSE PAD FITTING KIT

Further research has indicated the need to cater for a variety of nose shapes. Therefore, we have developed a range of different sized nose pads, Baby fit (B1, B2, B3), Asian fit (A1, A2, A3, A4, A5) and Wide fit (W1, W2, W3).


All nose pads are interchangeable on frames and the fitting kit stores all nine sizes of nose pad.

This kit will allow optometrists and opticians to custom fit the frames to each child, therefore, offering the best customer service.






THE RELATIONSHIP BETWEEN THE TYPES OF NOSE PADS AND FITTING RESULTS



If the frame sits too far away from a child's eyes, a wider nose pad should be used (e.g. replace a W1 for a W3 nose pad). If the frame sits too close to eyelashes then the above spacer kit can be used to create a greater space.


FOR INFANTS (B1,B2,B3), PLEASE REPLACE IT WITH THE LOWER NUMBER.

FUNCTION OF A4, A5 NOSE PADS



The A4 and A5 nose pads sets raise the eyeglass frame further up.

SPACER KIT



TOMATO GLASSES CASSETTE AND ACCESSORIES



- · Accessories are included in the Tomato Glasses cassette.
- \cdot The Baby frames have more accessories in their packs.
- · Very convenient for storing parts.

Contents

- \cdot frame \cdot strap \cdot extra nose pads \cdot extra tips
- \cdot extra tip screws \cdot driver \cdot fabric pouch \cdot user guide

EYE CATCHING PROMOTIONAL ITEMS

 $\ensuremath{\ensuremath{\%}}$ Contents and design are subject to change without any notice.



< 1set: 4ea >









OUR AWARDS

TOMATO GLASSES are an award winning specialized brand of children's glasses. We supply frames worldwide and we have a good reputation for design and quality from ophthalmologists and opticians around the world.



EYEWEAR OF TEAR 2016 AWARD(JAPAN)



KIDS DESIGN AWARD 2014



Optra Awards 2013 1st



Certificate of Good Design Selection



ODMA AWARDS 2017 (AUSTRALIA)



GOOD DESIGN (JAPAN)



AUSTRALIA O=MEGA AWARD



One million dollar, Prize for Top of export



Three million dollar, Prize for Top of export



Five million dollar, Prize for Top of export



ISO 9001 Certification



ISO 13485 Certification



The most Recommended Brand by Consumer



Toxin-Not Detected TEST REPORT(ear tips)



Toxin-Not Detected TEST REPORT(nose pads)



TEST REPORT(frame)
Toxin-Not Detected



Certificate Europe



No.10-1179157 Certificate of patent



No.10-1186259 Certificate of patent



Certificate of Trade Service Mark Registration



Cerificate Of Origin



No.30-1027921 Certificate of Design



No.30-0652187 Certificate of Design



No.30-0652189 Certificate of Design



Certificate of Trade Service Mark Registration



Certificate of Trademark Registration(JAPAN)



Registration(USA)



Certificate of Trademark Registration(INDIA)



Registration(Europe)



Certificate of Trademark Registration(China)



• Australia: TOMATO GLASSES AUSTRALIA (+61) 3-9448-8932

• Brazil : Vis-à-vis Importação (+55) 11-26797902

• Bulgaria : Milde Eood (+359) 888-811-508

• Canada: TOMATO GLASSES CANADA (+1) 778-954-7880

• Colombia : TIENDA OPTICA (+57) 909-485-5592

Czech Republic: TOMATO GLASSES Czech Republic (+420) 605747864
 HongKong: WING KI ENTERPRISE LIMITED (+852) 6014-6863
 Indonesia: PT.Optik Klinik Mata Nusantara (+62) 21-7581-6601
 Japan: Nagoya Spectacle Co., INC (+81) 52 262 6591

• Malaysia : TOMATO GLASSES (SINGAPORE) PTE LTD (+65) 9022-8960

Middle East: TOMATO GLASSES ME (+971) 4-252-4596
 New Zealand: TOMATO GLASSES NZ (+64) 21-869-300
 Philippines: TOMATO GLASSES PH (+63) 920-932-6578

• Singapore: TOMATO GLASSES (SINGAPORE) PTE LTD (+65) 9022-8960

• South Africa: IE Optics (PTY) LTD (+27) 11 485 0594

• Taiwan : ABBY Optical Limited (+886) 4 24711616 / GinDin international Co.,Ltd (+886) 2-2273-7785

Türkiye: MK OPTİK İÇ VE DIŞ TİCARET LİMİTED ŞİRKETİ (+90) 312-467-27-12
 United Kingdom: TOMATO GLASSES EU (+44) 1395-488188

• United States: TOMATO GLASSES USA (+1) 201-655-5758

• China: China



All employees at Tomato Glasses are committed to better vision for infants and children.

All of our employees are diligently working to deliver new ideas to improve vision for infants and children.

"Let every child see the world clearly."

"Every child should be able to study, to exercise, and to play freely without any restraints."

"Every child is precious, and we hope every child to grow up healthy." With these kinds of thoughts in mind, our employees are relentlessly striving together to make better products.

















HEAD OFFICE

29 Gudeok-ro 120 beon-gil, Seo-gu, Busan 49246, Korea
Tel: +82-51-244-8760, Fax: +82-51-244-8761
www.tomatoglasses.com

FACTORY

21, Haeyang-ro 57beon-gil, Youngdo-gu, Busan 49007, Republic of Korea

